

Module:
Techniques in Neuroscience

Week 2:
Electrophysiology: Looking at live neurons in action



Dr Jonathan Robbins

Topic 1:
An introduction to electrophysiology
Part 2 of 3

Part 2

Extracellular recording of electrical activity

1

Extracellular recording (ER):

- field potentials
- whole nerve activity
- multi-unit activity
- single unit activity
- multi-electrode arrays (MEAs)

2

Intracellular recording (IR):

- activity within single cells
- sharp electrodes
- patch suction electrodes

3

Single channel recording (SCR):

- recording activity of single ion channels
- patch clamp-type electrodes

Definitions to aid understanding

Definitions:



Field potential: this is the electric potential in the extracellular space around neurons.

Nerve: a bundle of axons.

Compound axon potential: the sum of the activity in a number of nerve fibers (or axons).

Click **Next** to continue

Extracellular recordings


There are five versions of extracellular recordings:

- field potentials
- whole nerve recordings
- multicellular (multi-unit) recordings
- single unit recordings
- MEAs

In each case

- the electrode is outside but close to the neurons
- the electrodes pick up only field potentials and low frequency filtered action potentials
- it is not possible to record V_m rest or post-synaptic potentials

Field potentials





O'keefe & Nadel (1978); The Scripps Research Institute (2008)

Whole nerve activity (1)

The compound action potential of a whole nerve can be recorded extracellularly.

Frog sciatic nerve action potential recorded with sodium as the permeant ion



Lilley & Robbins (1998), Robbins (unpublished data)

Whole nerve activity (2)


Galvani – a recreation of his experiment on frogs



The isolation of a whole nerve from a frog and recording of the action potential

These videos will comprise topic 2 for this week.

Whole nerve activity (3)



Different axons can be separated by the intensity of the stimulus and their conduction velocity.

Docherty et al. (2005)



Week 2 Electrophysiology: Looking at live neurons in action

Topic 1: An introduction to electrophysiology

9 of 16

Multi-unit activity

Multi-unit extracellular recordings can be done *in vivo*.



Robbins (unpublished data)

Week 2 Electrophysiology: Looking at live neurons in action


Topic 1: An introduction to electrophysiology

10 of 16

Single unit activity




An example of single unit recordings in the human medial temporal lobe.



Quiroga et al. (2009)

Multi-electrode arrays (1)

Multi-electrode array electrodes are bedded in the bottom of a dish.




Characteristics:

- 64 electrodes
- inert to cells
- pick up extracellular electrical activity

Neurons are grown inside these dishes in an incubator and monitor the electrical activity in a non-invasive manner.

Multi-electrode arrays (2)




Robbins (unpublished data)

Week 2 Electrophysiology: Looking at live neurons in action

Topic 1: An introduction to electrophysiology

13 of 16

Multi-electrode arrays (3)



Robbins (unpublished data)

Week 2 Electrophysiology: Looking at live neurons in action

Topic 1: An introduction to electrophysiology

14 of 16

Summary on extracellular recordings

Technique	Advantages	Disadvantages
Field potentials	Large number of neurons Record network activity <i>In vivo & in vitro</i>	Cannot identify individual cells Cannot monitor PSPs or Vm
Whole nerve	Large number of axons Can identify subgroups of axons by conduction velocity and stimulus threshold <i>In vivo & in vitro</i>	Cannot identify individual cells Cannot monitor PSPs or Vm
Multi-unit	Monitor a few cells simultaneously Characterise single axon or neuron <i>In vivo & in vitro</i>	Cannot monitor PSPs or Vm
Single units	Characterised single axon or neuron <i>In vivo & in vitro</i>	Difficult to do <i>in vivo</i> Cannot monitor PSPs or Vm
MEAs	Minimally invasive Large number of cells Record network activity Characterised single axon or neuron	Cannot target particular cells Not <i>in vivo</i> Cannot monitor PSPs or Vm

Click **Next** to continue

End of part 2