



Carl Hobbs

**Module:**  
**Techniques in Neuroscience**

Week 3:  
Immunohistochemistry: Preserving and studying cells of the brain


**Topic 1:**  
**An introduction to  
immunohistochemistry**  
Part 2 of 4

# Part 2

## Tissue processor (1)

**'Dip and dunk' machine:**  
12 stations

**How it works:**



Week 3 Immunohistochemistry: Preserving and studying cells of the brain


Topic 1: An introduction to immunohistochemistry

3 of 15

## Tissue processor (2)

**Enclosed tissue processor**

**How it works:**




Week 3 Immunohistochemistry: Preserving and studying cells of the brain

Topic 1: An introduction to immunohistochemistry


4 of 15

## Paraffin wax embedding station


### How it works:



## Embedding tissue in mould filled with paraffin wax



## Wax-embedded tissues: process




Week 3 Immunohistochemistry: Preserving and studying cells of the brain

Topic 1: An introduction to immunohistochemistry

7 of 15

## Sectioning




Week 3 Immunohistochemistry: Preserving and studying cells of the brain

Topic 1: An introduction to immunohistochemistry


8 of 15

## Microtomes: rotary (1)

## Benchtop rotary microtome




Sections are floated onto water in a bath ( $\sim 40^{\circ}\text{C}$ ).



Sections may then be separated to be mounted on slides individually or as ribbons.

## Microtomes: rotary (2)



## Microtomes: vibrating (1)


## Vibrating microtome (vibratome)

**Characteristics:**

- embedding media: agarose or gelatine
- section thickness: 50-500 micron thin sections


Sections are collected and are stained as free-floating and then mounted onto slides for microscopic examination.

## Microtomes: vibrating (2)



## Microtomes: freezing/cryostat

## Cryostat

**Characteristics:**

- used for sectioning frozen tissues
- temperature: -20°C
- cut sections can be immediately fixed (if frozen as fresh tissue) or stored
- vitreous water is hard enough for cutting sections
- embedding medium may be used if required

## Microtomes: sliding

## Sliding microtome

**Characteristics:**

- used to section frozen samples without the need for a relatively more expensive cryostat
- fitted to the benchtop
- tissue is kept frozen by blasting with CO<sub>2</sub> or solid CO<sub>2</sub>
- can produce 15-200 micron sections which are stained as free-floating sections

# End of part 2